

All masks are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Super soft ear loop makes comfortable wearing for long time

Nose Clip

Adjustable Nose Bridge,

Soft Elastic Ear Strap

High Elastic ear hook earring ears for a long time.

Water proof

Water proof test, more breathable

Fold Design

Easy to receive

Product Specifications:


Name	Disposable KF94 face mask
Fabric material	Layer 1: PP non woven (customized printing) Layer 2: high filtration paper Layer 3: high filtration paper Layer4: pp non woven
Style	4ply ,with round ear loop
Dimensions size	20 * 8 cm
Fabric weight	30gsm + 25gsm + 25gsm+25gsm, could be as your request
Color	Colorful mask
Certificate	CE, ISO13485, FDA
MOQ	20000 PCS
Packing design	can be as your design
Packing	50 pcs /box,2000 pcs /ctn, or as your requirement
Capability	1 million pcs per day
Delivery time	7-30 days
Applications	hospital, laboratory, clinic, food industry, dust-free workshop

The picture of KF94











Invisible Adjustable Soft Nose Clip
Better Fit Your Face
Barrier the Up Direction Air




High Strength Ear Loop
Comfortable Ware
No Hurt or Harm to Your Ear with
Long Period Wearing



Extra Soft Accessible
Skin Layer Inside soft
Non-Woven Protective Layer
No Harm to your Skin

HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.