

Product details

The property of medical facial masks

- 1) Polypropylene spun bond, soft and latex free
- 2) Adjustable aluminum or plastic nose piece, comfortable to wear
- 3) 3-Ply lace mask, high filtration
- 4) Elastic or tie on ear loop for your choice

How to wear disposable mouth cover

- 1) Flat mask
- 2) Push mask flat to the face
- 3) Hold the mask with left hand, and handle the lanyard around the ear portion with right hand
- 4) Hold the mask with right hand, and handle the lanyard around the ear portion with left hand.
- 5) Open the folds, and press the nose clip.

Packaging of surgical facemasks

Specification	Normal packing	Carton size	1*40 HQ container can hold (Pcs)
17.5*9.5cm, 3ply, 18+20+20G	50 pcs/box, 40 boxes/ctn	52*19*42cm	3,280,000 Pcs
17.5*9.5cm, 3ply, 25+25+25G	50 pcs/box, 40 boxes/ctn	51*37*30cm	2,400,000 Pcs
Remark: all can be special customized according to your actual requirements			

Product display



The regular shipment way



FAQ

Q: 1. What is your output capacity for your major range of the product?

A: Face mask capacity: 1000,000 pieces / day.

Coverall, Surgical gowns, lab coat, bed sheet capacity: 10,000 pieces/day.

Caps capacity: 200,000 pieces/day. Each 1*40HQ can be finished within 1 week.

Q:2. What is your MOQ of each products?

A: Face mask/caps/shoe covers/arm sleeves each MOQ: 50,000 pieces.

Clothes and sheets MOQ:5000 pieces