

□□□□□ □□□□□□□:

4layers **KF94 mask** with earloops,Adult 20*8CM,

The first layer: non-woven fabric (outer layer)

The second layer: melt blown cloth

The third layer: melt blown cloth

The fourth layer: non-woven fabric (inner layer)

□□□□ □□□□

Adjustable Nose Bridge,

Soft Elastic Ear Strap


High Elastic ear hook earring ears for a long time.



□□□□ □□

Water proof test, more breathable

Fold Design

Easy to receive




	Invisible Adjustable Soft Nose Clip Better Fit Your Face Barrier the Up Direction Air		High Strength Ear Loop Comfortable Ware No Hurt or Harm to Your Ear with Long Period Wearing		Extra Soft Accesible Skin Layer Inside soft Non-Woven Protective Layer No Harm to your Skin
---	---	---	---	---	--

□□□□□ □□□

□□□	Disposable KF94 face mask
□□□□	Nonwoven face mask

□□□□	Black,Blue , Red, Green, Yellow
□□□□□ □□□□	Hubei, China (Mainland)
□□□□	1.4-ply protection,highly breathable
	2.PFE>95%,BFE>95%
	3.Style: Adjustable nose bridge Tie on
	4.high filtration efficiency
Ply Type	4 ply
□□□□	earloop
Applications	hospital,laboratory,food industry,dust free workshop,electronic manufacturers and so on
□□□□□□	20pcs/box,2000pcs/ctn,or as your requirements
Remark	all can be special customized according to your actual

HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.



