



3D	3D
BFE PFE.	bfe≥95%
Material	PP + Nonwoven fabric + Nonwoven fabric + PP
Color	Colorful / Customized
Size	50G + 25G + 25G + 30G
Dimensions	20x8cm
Quantity	50 40 / bags / 2000 pcs
Lead time	7-15 days
Payment	t / t 1 / c / x-transfer / paypal
Service	OEM / ODM
Standard	GB2626-2019.

orient Honest Medical International is a professional manufacturer of various types of masks.

* 4 types of masks available

orient Honest Medical International is a professional manufacturer of various types of masks *

orient Honest Medical International is a professional manufacturer of various types of masks *

orient Honest Medical International *

orient Honest Medical International is a professional manufacturer of various types of masks *



□□□□□□ □□□□□□ □□□□□□ □□□□□□ □□□□□ □□□□ □□ □□□□

HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.



Company Profile



Factory outside



Material warehouse



Dust free mask workshop



Dust free sewing workshop



Dust free packing workshop



Finished products warehouse



CE



ISO 13486



FDA