

All masks are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Super soft ear loop makes comfortable wearing for long time

Nasenklammer

Einstellbare Nasenbrücke,

Weiche elastische Ohrgurt

Hohe elastische Ohrhaken-Ohren-Ohren lange Zeit.

Wasserdicht

Water proof test, more breathable

Faltendesign


Leicht empfangen

Produktspezifikationen:

Name	Einweg-KF94-Gesichtsmaske
Fabric material	Layer 1: PP non woven (customized printing) Layer 2: high filtration paper Layer 3: high filtration paper Layer4: pp non woven
Stil	4ply ,with round ear loop
Dimensions size	20 * 8 cm
Fabric weight	30gsm + 25gsm + 25gsm+25gsm, could be as your request
Farbe	Colorful mask
Zertifikat	CE, ISO13485, FDA
Moq.	20000 PCS
Packing design	can be as your design
Verpackung	50 pcs /box,2000 pcs /ctn, or as your requirement
Fähigkeit	1 million pcs per day
Lieferzeit	7-30 days
Anwendungen	hospital, laboratory, clinic, food industry, dust-free workshop

Proben können wie gewünscht bereitgestellt werden. Willkommen Ihre Anfrage.

The picture of KF94



HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.





