

## Product Description

Name	Einweg-KF94-Gesichtsmaske
Material	Vlies-Gesichtsmaske
Farbe	Schwarz, Blau, Rot, Grün, Gelb
Herkunftsort	Hubei, China (Festland)
Spezifikationen	1.4-PLY-Schutz, hoch atmungsaktiv 2.PFE> 95%, BFE> 95% 3.style: Einstellbare Nasenbrücke Bindung an 4. Hohe Filtrationseffizienz
PLY-Typ	4 Ply
Stil	oldopf
Anwendungen	Krankenhaus, Labor, Lebensmittelindustrie, staubfreier Werkstatt, elektronisch Hersteller und so weiter
Verpackung	20pcs / Box, 2000pcs / ctn oder als Ihre Anforderungen
Hinweis	1. Dieses Produkt wird von Ethylenoxid sterilisiert 2.Dieselbar
Bemerkung	Alle können nach Ihrem tatsächlichen Spezial sein

## HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

## Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.





