



Nombre	Mascarilla de cara desechable de los niños.
Material	Tela no tejida
Color	impresión
Lugar de origen	Hubei, China (Mainland)
Especificaciones	Protección de 1.4-capa, altamente transpirable. 2.PFE> 95%, BFE> 95% 3.Style: Lazo de puente de nariz ajustable en 4. alta eficiencia de filtración
Tipo de capa	4 capitán
Estilo	oreja
Aplicaciones	Hospital, laboratorio, industria alimentaria, Taller de polvo, electrónica. Fabricantes y así sucesivamente
Embalaje	20pcs / caja, 2000pcs / ctn, o como sus requisitos
Nota	1. Este producto se esteriliza con oxid de etileno. 2.Disposición
Observación	Todos pueden ser personalizados especiales de acuerdo a su real.

Las muestras se pueden proporcionar según lo solicitado. Bienvenido a su consulta.



HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.







