


Description du produit

Nom	Masque de visage KF94 jetable
Matériau	masque de visage non tissé
couleur	Noir, Bleu, Rouge, Vert, Jaune
Lieu d'origine	Hubei, Chine (continentale)
Spécifications	Protection de 1,4 plis, très respirante 2.PFE> 95%, BFE> 95% 3.Style: cravate de pont de nez ajustable sur 4.High Filtration Efficacité
Type de pli	4 pli
style	poursuivre
Applications	Hôpital, laboratoire, industrie alimentaire, atelier sans poussière, électronique Les fabricants et ainsi de suite
Emballage	20pcs / boîte, 2000PCS / CTN ou selon vos exigences
Note	1.Ce produit est stérilisé par oxyde d'éthylène 2.Disposable
remarque	Tout peut être spécial personnalisé en fonction de votre réel

HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.





