




| | |
|-----------------------|--|
| Nama Produk | Masker Wajah KF94 sekali pakai |
| BFE. | BFE≥95% |
| Bahan | PP + Melt Blown + Melt Blown + PP |
| Warna | Berlipat Warna (Dapat Disesuaikan) |
| Gaya | Loop telinga |
| Bobot | 50g + 25g + 25g + 30g, atau sebagai permintaan Anda |
| Ukuran | 20x8cm untuk orang dewasa, kontak untuk ukuran anak-anak |
| Sedang mengemas | 50 pcs / bag, 40 bags / karton, 2000 pcs / karton |
| Waktu pengiriman | 7-15 hari |
| Syarat pembayaran | T / T, L / C, Paypal, dll |
| Melayani | OEM, sampel disediakan, penyediaan spesifikasi produk |
| Pengujian sertifikasi | GB2626-2019. |

Sampel dapat disediakan sesuai permintaan. Selamat datang pertanyaan Anda.



HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.