

Spesifikasi Produk:

4Layers KF94 Mask dengan earloops, dewasa 20*8cm,

Lapisan pertama: kain non-anyaman (lapisan luar)

Lapisan kedua: melelehkan kain yang ditiupkan

Lapisan ketiga: melelehkan kain yang melebur

Lapisan keempat: kain non-anyaman (lapisan dalam)

Klip hidung

Jembatan hidung yang bisa disesuaikan,

Tali telinga elastis lembut

Telinga telinga telinga elastis tinggi untuk waktu yang lama.




Deskripsi Produk

Nama	Masker wajah Kf94 sekali pakai
Bahan	Topeng Wajah Nonwoven
Warna	Hitam, biru, merah, hijau, kuning
Tempat asal	Hubei, Cina (Daratan)
Spesifikasi	1. 4-ply Protection, sangat bernapas
	2. pfe > 95%, BFE > 95%
	3. Gaya: Dasi Jembatan Hidung yang Dapat Disesuaikan
	4. Efisiensi filtrasi tinggi

Jenis lapis	4 ply
Gaya	earloop
Aplikasi	rumah sakit, laboratorium, industri makanan, lokakarya bebas debu, elektronik produsen dan sebagainya
Sedang mengemas	20pcs/box, 2000pcs/ctn, atau sebagai kebutuhan Anda
Komentar	Semua dapat disesuaikan khusus sesuai dengan aktual Anda

HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.



