



Deskripsi Produk:

Nama	Anak -anak sekali pakai KF94 masker
Bahan	Non-anyaman masker
Warna	Pola warna
Tempat asal	Hubei, Cina (Daratan)
Spesifikasi	1. Sangat bernafas 2. PFE > 95%, BFE > 95% 3. Gaya: Jembatan Hidung yang Dapat Disesuaikan, Ear_Loop 4. 4-ply Perlindungan, efisiensi filtrasi yang tinggi
Jenis lapis	4 ply
Gaya	Ear_loop
Aplikasi	rumah sakit, laboratorium, industri makanan, lokakarya bebas debu, elektronik produsen dan sebagainya
Sedang mengemas	20pcs/box, 2000pcs/ctn, atau sebagai kebutuhan Anda
Catatan	1. Produk ini disterilkan dengan etilena oksid 2. sekali pakai
Komentar	Semua dapat disesuaikan khusus sesuai dengan aktual Anda

Sampel dapat disediakan seperti yang diminta. Selamat datang pertanyaan Anda.



HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.



