


Nama Wajah KF94 sekali pakai [masker](#)
Material Masker wajah nonwoven
Warna Hitam, biru, merah, hijau, kuning
tempat asal Hubei, China (Mainland)
Spesifikasi 1. [4 lapis](#) perlindungan, sangat bernapas
2. $Pfe > 95\%$, $BFE > 95\%$
3. Style: ikat jembatan hidung yang dapat disesuaikan
4. efisiensi filtrasi tinggi
jenis ply 4 lapis
Gaya earloop.
Aplikasi Rumah Sakit, Laboratorium, Industri Makanan, Lokakarya Bebas Debu, Elektronik produsen dan sebagainya
Packing 20 pcs / box, 2000pcs / ctn, atau sebagai kebutuhan Anda
Catatan 1. Produk ini disterilkan oleh ethylene oxid
2. Aposisi
komentar Semua bisa disesuaikan khusus sesuai dengan yang sebenarnya

HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.





