



nome del prodotto	Maschera viso kf94 monouso
BFE.	BFE≥95%
Materiale	PP + Melt Blown + Melt Blown + PP
Colore	Colore classificato (personalizzabile)
Stile	Ear-Loop.
Peso	50g + 25g + 25g + 30g, o come vostra richiesta
Dimensione	20x8 cm per adulti, contatto per le dimensioni dei bambini
Imballaggio	50pcs / bag, 40Bags / cartone, 2000pcs / cartone
Tempo di consegna	7-15 giorni
Modalità di pagamento	T / T, L / C, PayPal, ecc
Servizio	OEM, campioni forniti, fornitura di specifiche del prodotto
Certificazione di test	GB2626-2019.

I campioni possono essere forniti come richiesto. Benvenuto la tua richiesta.



# HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

## Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.