
nome del prodotto	Maschera viso kf94 monouso
BFE.	BFE≥95%
Materiale	PP + Melt Blown + Melt Blown + PP
Colore	Colore classificato (personalizzabile)
Stile	Ear-Loop.
Peso	50g + 25g + 25g + 30g, o come vostra richiesta
Dimensione	20x8 cm per adulti, contatto per le dimensioni dei bambini
Imballaggio	50pcs / bag, 40Bags / cartone, 2000pcs / cartone
Tempo di consegna	7-15 giorni
Modalità di pagamento	T / T, L / C, PayPal, ecc
Servizio	OEM, campioni forniti, fornitura di specifiche del prodotto
Certificazione di test	GB2626-2019.

I campioni possono essere forniti come richiesto. Benvenuto la tua richiesta.



## **HOW TO USE**



- Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

## Precautions for Use

- 1. Make sure not to place towel, tissue, etc. inside of the mask.
- 2. Do not use if the inside of the mask is contaminated.
- 3. Do not touch the surface of the mask after wearing.
- 4. This mask is disposable. Do not wash your mask.