



Nome	Maschera per il viso KF94 monouso per bambini
Materiale	Tessuto non tessuto
Colore	stampa
Luogo d'origine	Hubei, China (Mainland)
Specifiche	Protezione da 1.4-Ply, altamente traspirante
	2.pfe> 95%, BFE> 95%
	3.Style: cravatta regolabile del ponte del naso
	4. Efficienza di filtrazione
Type.	4 strati
Stile	in primo piano
Applicazioni	Ospedale, Laboratorio, Industria alimentare, Laboratorio privo di polvere, Elettronica produttori e così via
Imballaggio	20pcs / box, 2000pcs / ctn o come vostre esigenze
Nota	1. Questo prodotto è sterilizzato da etilene ossido
	2.Disponibile
Osservazione	Tutti possono essere speciali personalizzati in base al tuo attuale

I campioni possono essere forniti come richiesto. Benvenuto la tua richiesta.



HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.







