

## [4ply nonwoven medical kf94 face mask for adults](#)

**All masks are made of the highest quality materials, providing excellent protection and superior comfort:**

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Super soft earloop makes comfortable wearing for long time

## [4ply nonwoven medical kf94 face mask for adults](#)

|                        |  |
|------------------------|--|
| <b>Fabric Material</b> | Layer 1: PP non woven (customized printing)<br>Layer 2: high filtration paper<br>Layer 3:high filtration paper<br>Layer4: pp non woven |
| <b>Style</b>           | 4ply ,with round ear loop  |
| <b>Dimensions size</b> | 20*8cm   |
| <b>Fabric weight</b>   | 30gsm + 25gsm + 25gsm+25gsm, could be as your request  |
| <b>Color</b>           | Colorful mask  |
| <b>Certificate</b>     | CE, ISO13485, FDA  |
| <b>MOQ:</b>            | 20000pcs   |
| <b>Packing design</b>  | can be as your design  |
| <b>Packing</b>         | 50 pcs/box,2000 pcs/ctn, or as your requirement  |
| <b>Capability</b>      | 1 million pcs per day  |
| <b>Delivery time</b>   | 7-30 days  |
| <b>Applications</b>    | hospital, laboratory, clinic, food industry, dust-free workshop  |

## [The picture of 4ply nonwoven medical kf94 face mask for adults](#)









# HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

## Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.