

4ply nonwoven medical printing kf94 face mask for adults


All masks are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Super soft earloop makes comfortable wearing for long time


4ply nonwoven medical printing kf94 face mask for adults

| | |
|-------------------|---|
| Fabric Material | Layer 1: PP non woven (customized printing) Layer 2: high filtration paper Layer 3: high filtration paper Layer4: pp non woven |
| Stile | 4ply ,with 3.2mm flat ear loop |
| Dimensions size | 20*8cm |
| Fabric weight | 30gsm + 25gsm + 25gsm+25gsm, could be as your request |
| Colore | printing pattern mask |
| Certificato | CE, ISO13485, FDA |
| MOQ: | 20000pcs |
| Packing design | can be as your design |
| Imballaggio | 50 pcs/box, 2000 pcs/ctn, or as your requirement |
| Capacità | 1 million pcs per day |
| Tempo di consegna | 7-30 days |
| Applicazioni | hospital, laboratory, clinic, food industry, dust-free workshop |

The picture of 4ply nonwoven medical printing kf94 face mask for adults



HOW TO USE



- 1** Remove mask from packaging.
- 2** Wear the mask to cover your nose and chin.
- 3** Hook each strap over your ear to lock them in place
- 4** Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.





