



型号	KN94
b	BFE95
材料	PP + 无纺布 + 无纺布 + PP.
重量	50g + 25g + 25g + 30g
尺寸	20x8cm
数量	50个/包 40个/包 2000pcs /包
保质期	7-15天
支付方式	T / T / L / C / PayPal
其他	OEM
标准	GB2626-2019

KN94



	<p>Invisible Adjustable Soft Nose Clip Better Fit Your Face Barrier the Up Direction Air</p>		<p>High Strength Ear Loop Comfortable Wear No Hurt or Harm to Your Ear with Long Period Wearing</p>		<p>Extra Soft Accessible Skin Layer inside soft Non-Woven Protective Layer No Harm to your Skin</p>
--	--	--	---	--	---

# HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

## Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.