


|            |                                                                                                        |
|------------|--------------------------------------------------------------------------------------------------------|
| Model      | Model KF94                                                                                             |
| Material   | Non-woven                                                                                              |
| Color      | White                                                                                                  |
| Size       | One size                                                                                               |
| Features   | 1.4-layer non-woven fabric<br>2.PFE > 95% BFE > 95%<br>3. Ultra-soft non-woven fabric<br>4. Ultra-soft |
| Quantity   | 4                                                                                                      |
| Weight     | 1.5g                                                                                                   |
| Usage      | For daily use, travel, and outdoor activities                                                          |
| Package    | 20 PCS / BOX 2000PCS / CTN                                                                             |
| Storage    | Store in a cool, dry place                                                                             |
| Shelf life | 12 months                                                                                              |

1. Ultra-soft non-woven fabric




Invisible Adjustable Soft Nose Clip  
Better Fit Your Face  
Barrier the Up Direction Air

High Strength Ear Loop  
Comfortable Wear  
No Hurt or Harm to Your Ear with Long Period Wearing

Extra Soft Accessible Skin Layer inside soft Non-Woven Protective Layer  
No Harm to your Skin

# HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

## Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.













