

**All mask. It is made of high quality material and offers excellent protection and excellent comfort.**

-Son-Bright Materials are useful for excluding smooth and skin stimulation 95%

**Product Specifications:**

**4layers KF94 mask with earloops,20\*8cm (Accept customized size)**

The first layer: SS grade 50g non-woven fabric (outer layer)

The second layer: 99 grade 25 grams of melt blown cloth

The third layer: 99 grade 25 grams of melt blown cloth

The fourth layer: SS grade 30g non-woven fabric (inner layer)

Nose Clip

Adjustable Nose Bridge,

Soft Elastic Ear Strap


High Elastic ear hook earing ears for a long time.

Water proof

Water proof test, more breathable


Fold Design

Easy to receive



	<p>Invisible Adjustable Soft Nose Clip Better Fit Your Face Barrier the Up Direction Air</p>		<p>High Strength Ear Loop Comfortable Ware No Hurt or Harm to Your Ear with Long Period Wearing</p>		<p>Extra Soft Accessible Skin Layer Inside soft Non-Woven Protective Layer No Harm to your Skin</p>
---	--	---	---	---	---

# HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

## Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.









orient Honest  
MEDICAL INTERNATIONAL

