



Product Name	Disposable KF94 face mask
BFE	BFE≥95%
Material	PP+ Melt Blown+Melt Blown+ PP
Color	color-graded (Customizable)
Style	Ear-loop
Weight	50g+25g+25g+30g, or as your request
Size	20x8cm for Adults, contact for kids' size
Packing	50pcs/bag, 40bags/carton, 2000pcs/carton
Delivery Time	7-15 days
Terms of Payment	T/T, L/C, PayPal, etc
Service	OEM, Samples provided, Provision of Product Specification
Testing Certification	GB2626-2019

Samples can be provided as requested. Welcome your enquiry.



HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.