

All masks are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Super soft ear loop makes comfortable wearing for long time

Nose Clip

Adjustable Nose Bridge,

Soft Elastic Ear Strap

High Elastic ear hook earing ears for a long time.

Water proof

Water proof test, more breathable

Fold Design


Easy to receive

Product Specifications:

| | |
|-----------------|---|
| Name | Disposable KF94 face mask |
| Fabric material | Layer 1: PP non woven (customized printing) Layer 2: high filtration paper Layer 3: high filtration paper Layer4: pp non woven |
| Style | 4ply ,with round ear loop |
| Dimensions size | 20 * 8 cm |
| Fabric weight | 30gsm + 25gsm + 25gsm+25gsm, could be as your request |
| Color | Colorful mask |
| Certificate | CE, ISO13485, FDA |
| MOQ | 20000 PCS |
| Packing design | can be as your design |
| Packing | 50 pcs /box,2000 pcs /ctn, or as your requirement |
| Capability | 1 million pcs per day |
| Delivery time | 7-30 days |
| Applications | hospital, laboratory, clinic, food industry, dust-free workshop |

Samples can be provided as requested. Welcome your enquiry.

The picture of KF94



HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.





