

## Product details

### The indication of isolation mouth mask

- 1) The mask is a fluid resistant and single-use mask
- 2) Protect others from bacteria exhaled by the wearer
- 3) Reduces exposure of the wearer to splashes of blood and bodily fluids.
- 4) It's light weight, and designed to provide comfort and easy breathing

### How to use wind protection face mask

- 1) Place mask over nose and mouth, and place an ear loop around each ear
- 2) Pull bottom of mask under chin and mold the nose piece around the nose
- 3) The mask should cover face from the nose to underneath the chin
- 4) Discard mask after use, do not reuse masks

### Packaging of face mask surgical disposable 3Ply


Specification	Normal packing	Carton size	1*40 HQ container can hold (Pcs)
17.5*9.5cm,3ply, 18+20+20G	50 pcs/box,40 boxes/ctn	52*19*42cm	3,280,000 Pcs
17.5*9.5cm,3ply, 25+25+25G	50 pcs/box,40 boxes/ctn	51*37*30cm	2,400,000 Pcs

Remark: all can be special customized according to your actual requirements

## Product display



## Our business process



## The regular shipment way

