

CE ISO approved China face mask higienic products medical supplies manufactuer

Property of face mask

- 1) High filtration micro-fiber filter media with a BFE rating 99%.
- 2) Each mask features three layers, with the innermost layer soft and smooth for maximum comfort.
- 3) Elastic earloop style, fluid resistant and latex-free. Help protect against the Swine flu and other epidemics.
- 4) Special size for the ladies: 145×90 MM

Packaging of higienic products face mask

Specification	Normal packing	Carton size	1*40 HQ container can hold (Pcs)
3-ply, 25+25+25G	50 pcs/box,40 boxes/ctn	52*31*34cm	2,440,000 Pcs
Remark: all can be special customized according to your actual requirements			

How to wear the medical supplies face mask

- 1) Flat mask
- 2) Push mask flat to the face
- 3) Hold the mask with left hand, and handle the lanyard around the ear portion with right hand
- 4) Hold the mask with right hand, and handle the lanyard around the ear portion with left hand.
- 5) Open the folds, and press the nose clip.

Pictures of face mask



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

