

## **China exporters for industries mouth mask , simple ear loop face mask , fashionable surgical masks**

### **Properties of mouth mask**

- 1) High filtration:3-ply pleated design offers,0.05 micron filtration at greater than 99%.
- 2) Slim aluminum wire fits comfortably over the nose,and maintain maximum filtration at the same time.
- 3) Tie and ear loop style: two pairs of polypropylene ties insure a comfortable fit; ear loop style attaches fast and with elastic loops.

### **Specifications of ear loop face mask**

Material	PP non woven + high filtration paper + PP non woven
Fabric weight	Could be as your request
Ear loop	Tie on or ear loop
Size	17.5*9.5cm, or as your requirements
Color	white/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black
Packing	50 pcs/box,2000 pcs/ctn, or as your requirement
Packing design	all inner box and outer carton printing can be as your design
Applications	hospital,laboratory,food industry,dust-free workshop,electronic manufactures and so on

### **Storage of ear loop face mask**

- 1) The product will maintain its properties when it is stored in a dry place at temperatures between 10 ° and 30 ° C.
- 2) Not exposed to direct sunlight.

### **How to wear our fashionable surgical masks**

- 1) Flat mask
- 2) Push mask flat to the face
- 3) Hold the mask with left hand,and handle the lanyard around the ear portion with right hand
- 4) Hold the mask with right hand,and handle the lanyard around the ear portion with left hand.
- 5) Open the folds,and press the nose clip.

### **Picture of mouth mask**

