
#### 3ply non woven disposable printed face masks for adults

# All <u>masks</u> are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Super soft earloop makes comfortable wearing for long time





## **Product Description**

Name	Disposable <u>KF94 face mask</u>	
Material	Nonwoven face mask	
Color	Black,Blue , Red, Green, Yellow	
Place of Origin	Hubei, China (Mainland)	
Specifications	1.4-ply protection, highly breathable	
	2.PFE>95%,BFE>95%	
	3.Style: Adjustable nose bridge Tie on	
	4.high filtration efficiency	
Ply Type	4 ply	
Style	earloop	
Applications	hospital,laboratory,food industry,dust free workshop,electronic	
	manufacturers and so on	
Packing	20pcs/box,2000pcs/ctn,or as your requirements	

Note	1. This product is sterilized by ethylene oxid	
	2.Disposable	
Remark	all can be special customized according to your actual	

#### **Product Specifications:**

**4layers KF94 mask with earloops,**20\*8cm (Accept customized size)

The first layer: SS grade 50g non-woven fabric (outer layer) The second layer: 99 grade 25 grams of melt blown cloth The third layer: 99 grade 25 grams of melt blown cloth

The fourth layer: SS grade 30g non-woven fabric (inner layer)

### **Nose Clip**

Adjustable Nose Bridge,

### **Soft Elastic Ear Strap**

High Elastic ear hook earing ears for a long time.


#### Water proof

Water proof test, more breathable

#### **Fold Design**

Easy to receive

# **HOW TO USE**



- Remove mask from packaging.
- Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

## **Precautions for Use**

- 1. Make sure not to place towel, tissue, etc. inside of the mask.
- 2. Do not use if the inside of the mask is contaminated.
- 3. Do not touch the surface of the mask after wearing.
- 4. This mask is disposable. Do not wash your mask.









