

## **All masks are made of the highest quality materials, providing excellent protection and superior comfort:**

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Super soft ear loop makes comfortable wearing for long time

### **Nose Clip**

Adjustable Nose Bridge,

### **Soft Elastic Ear Strap**

High Elastic ear hook earring ears for a long time.

### **Water proof**

Water proof test, more breathable

### **Fold Design**

Easy to receive

### **Product Specifications:**

Name	Disposable KF94 face mask
<b>Fabric material</b>	Layer 1: PP non woven (customized printing) Layer 2: high filtration paper Layer 3: high filtration paper Layer4: pp non woven
<b>Style</b>	4ply ,with round ear loop
<b>Dimensions size</b>	20 * 8 cm
Fabric weight	30gsm + 25gsm + 25gsm+25gsm, could be as your request
Color	Colorful mask
Certificate	CE, ISO13485, FDA
MOQ	20000 PCS
Packing design	can be as your design
Packing	50 pcs /box,2000 pcs /ctn, or as your requirement
Capability	1 million pcs per day
Delivery time	7-30 days
Applications	hospital, laboratory, clinic, food industry, dust-free workshop

### **The picture of KF94**









Invisible Adjustable Soft Nose Clip  
Better Fit Your Face  
Barrier the Up Direction Air




High Strength Ear Loop  
Comfortable Wear  
No Hurt or Harm to Your Ear with  
Long Period Wearing




Extra Soft Accessible  
Skin Layer Inside soft  
Non-Woven Protective Layer  
No Harm to your Skin





# HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

## Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.