

## **Disposable face mask with earloop , disposable hospital mask , disposable masks factory**

### **Quick details for disposable hospital mask**

- 1) Layer: 3ply, PP non woven + high filtration paper + PP non woven
- 2) Size: 19.5X7.5cm
- 3) Color: white/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black, we have 11 different colors totally.
- 4) Ear loop: flat elastic, round elastic, or tie on ear loop
- 5) Nose bar: Plastic or metal

### **Our disposable masks can achieve following quality standard**

- 1) Conform to EN14683: 2005
- 2) Bacterial filtration efficiency (BFE%) > 99%.
- 3) TYPE IIR standard, resistance to penetration by synthetic blood-120mmHg.

### **How to wear disposable face mask with earloop**

- 1) Flat the disposable masks
- 2) Push disposable masks flat to the face
- 3) Hold the mask with left hand, and handle the lanyard around the ear portion with right hand
- 4) Hold the mask with right hand, and handle the lanyard around the ear portion with left hand.
- 5) Open the folds, and press the nose clip.

### **Packaging of disposable hospital mask**

Specification	Normal packing	Carton size	1*40 HQ container can hold (Pcs)
17.5*9.5cm, 3ply, 18+20+20G	50 pcs/box, 40 boxes/ctn	52*19*42cm	3,280,000 Pcs
17.5*9.5cm, 3ply, 25+25+25G	50 pcs/box, 40 boxes/ctn	51*37*30cm	2,400,000 Pcs
Remark: all can be special customized according to your actual requirements			

### **Other information of disposable face mask with earloop**

MOQ	100,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Supply ability	1,000,000 pcs/day, 1*40 HQ container can be finished in 15 days
Delivery time	Normally 30 days after confirmed the order
Sample	Sample is free, can be prepared well in 2 days

### **Pictures of disposable face mask with earloop**