


Nome do Produto	Máscara de rosto KF94 descartável
Bfe.	BFÉ≥95%
Material	Pp + derretimento soprado + derrete soprado + pp
Cor	Cor-classificado (personalizável)
Estilo	Loop de orelha
Peso	50g + 25g + 25g + 30G, ou como seu pedido
Tamanho	20x8cm para adultos, contato para o tamanho das crianças
Embalagem	50 pcs / saco, 40bags / caixa, 2000 pcs / caixa
Prazo de entrega	7-15 dias
Termos de pagamento	T / t, l / c, paypal, etc
Serviço	OEM, amostras fornecidas, fornecimento de especificação do produto
Certificação de teste	Gb2626-2019.

As amostras podem ser fornecidas conforme solicitado. Bem-vindo sua pergunta.



HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.