

Descrição do produto

nome

Máscara facial KF94 descartável

Material

máscara facial não tecida

cor

preto, azul, vermelho, verde, amarelo

Local de origem

Hubei, China (Mainland)

Especificações

Proteção de 1.4-PLY, altamente respirável

2.PFE> 95%, BFE> 95%

3.Style: laço de ponte de nariz ajustável em

4.High eficiência de filtração

tipo de ply

4 camadas

estilo

lagoa

Aplicações

Hospital, laboratório, indústria alimentícia, oficina livre de poeira, eletrônica
fabricantes e assim por diante

embalagem

20 pcs / caixa, 2000 pcs / ctn, ou como suas necessidades

nota


1.Este produto é esterilizado pela oxid do etileno

2.Disponível

Observação

Todos podem ser especiais personalizados de acordo com o seu

HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.





