

#### 4ply nonwoven medical printing kf94 face mask for adults


All masks are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Super soft earloop makes comfortable wearing for long time


#### 4ply nonwoven medical printing kf94 face mask for adults

Fabric Material	Layer 1: PP non woven (customized printing) Layer 2: high filtration paper Layer 3: high filtration paper Layer4: pp non woven
Estilo	4ply ,with 3.2mm flat ear loop
Dimensions size	20*8cm
Fabric weight	30gsm + 25gsm + 25gsm+25gsm, could be as your request
Cor	printing pattern mask
Certificado	CE, ISO13485, FDA
MOQ:	20000pcs
Packing design	can be as your design
Embalagem	50 pcs/box, 2000 pcs/ctn, or as your requirement
Capacidade	1 million pcs per day
Tempo de entrega	7-30 days
Formulários	hospital, laboratory, clinic, food industry, dust-free workshop

The picture of 4ply nonwoven medical printing kf94 face mask for adults



## HOW TO USE



- 1** Remove mask from packaging.
- 2** Wear the mask to cover your nose and chin.
- 3** Hook each strap over your ear to lock them in place
- 4** Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

## Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.





