

Todas as máscaras são feitas de materiais da mais alta qualidade, proporcionando excelente proteção e conforto superior:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Super soft ear loop makes comfortable wearing for long time

Nose Clip

Adjustable Nose Bridge,

Soft Elastic Ear Strap

High Elastic ear hook earring ears for a long time.

Water proof

Water proof test, more breathable

Fold Design

Easy to receive

Product Specifications:

Name	Disposable KF94 face mask
Fabric material	Layer 1: PP non woven (customized printing) Camada 2: Papel de filtração alta Layer 3: high filtration paper Layer4: pp non woven

Estilo

>>
4ply, com loop redondo de ouvido

Tamanho das dimensões

>>
20 * 8 cm

Peso do tecido

>>
30gsm + 25gsm + 25gsm + 25gsm, pode ser como sua solicitação

cor

>>
Máscara colorida

Certificado

>>
CE, ISO13485, FDA

MOQ

>>
20000 PCs

Design de embalagem

>>
pode ser como seu design

embalagem

>>
50 PCs /Box, 2000 PCs /CTN, ou como seu requisito

Capacidade

>>
1 milhão de PCs por dia

Tempo de entrega

>>
7-30 dias

Aplicativos

>>
Hospital, laboratório, clínica, indústria de alimentos, oficina sem poeira

>

A imagem do KF94>











Invisible Adjustable Soft Nose Clip
Better Fit Your Face
Barrier the Up Direction Air



High Strength Ear Loop
Comfortable Ware
No Hurt or Harm to Your Ear with
Long Period Wearing



Extra Soft Accessible
Skin Layer Inside soft
Non-Woven Protective Layer
No Harm to your Skin

HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.