

4ply nonply medical printing kf94 mask


Technical specifications for 4ply nonply medical printing kf94 mask:

- Material specifications
- Performance requirements
- Super Soft Earloop

4ply nonply medical printing kf94 mask

Material	Layer 1: PP (non-woven) Layer 2: Non-woven fabric Layer 3: Non-woven fabric Layer 4: PP (non-woven)
Structure	4ply non-woven fabric 3.2mm
Dimensions	20*8 cm
Weight	30GSM + 25GSM + 25GSM + 25GSM
Standards	CE, ISO13485, FDA
Moq	20000pcs
Lead time	Standard lead time
Packing	50 pieces/bag, 2000 pcs/ctn
Color	White
Shelf life	7-30 days
Storage	Store in a cool, dry place, away from sunlight and moisture.

4ply Nonply KF94 FACE MASK



HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.





