


Product Description:


Material	Non-woven fabric with <b>KF94</b> filter
Color	Black
Size	Standard
Origin	Hubei, China (Guangdong)
Features	<ol style="list-style-type: none"> <li>High filtration efficiency</li> <li>PFE &gt; 95%, BFE &gt; 95%</li> <li>Adjustable: Includes <b>EAR_LOOP</b></li> <li><b>4-Ply</b> non-woven fabric</li> </ol>
Ply	4 Ply
Attachment	EAR_LOOP
Usage	Indoor and outdoor use, suitable for daily protection, travel, and public places.
Quantity	20pcs/box, 2000pcs/ctn
Notes	<ol style="list-style-type: none"> <li>Check the fit and ensure it covers the nose and mouth.</li> <li>Do not touch the front of the mask.</li> </ol>
Storage	Store in a clean, dry place.

Product Features and Benefits:



<p>Invisible Adjustable Soft Nose Clip Better Fit Your Face Barrier the Up Direction Air</p>	<p>High Strength Ear Loop Comfortable Ware No Hurt or Harm to Your Ear with Long Period Wearing</p>	<p>Extra Soft Accessible Skin Layer inside soft Non-Woven Protective Layer No Harm to your Skin</p>
--	---	---

# HOW TO USE



- 1 Remove mask from packaging.
- 2 Wear the mask to cover your nose and chin.
- 3 Hook each strap over your ear to lock them in place
- 4 Press the nose piece with the fingers of both hands so that the mask is in close contact with your face.

## Precautions for Use

1. Make sure not to place towel, tissue, etc. inside of the mask.
2. Do not use if the inside of the mask is contaminated.
3. Do not touch the surface of the mask after wearing.
4. This mask is disposable. Do not wash your mask.



