

## **Xiantao supplier for disposable PP heavy weight patient slide sheet CE approved**

### **The specification for heavy weight patient slide sheet**

Material	polypropylene
Color	blue
Style	Ultrasonic heat sealing ,with hole in the length side
Size	195cm*85cm (can be made according to client's demand)
Application	Hospital, emergency center, clinic
MOQ	5,000 Pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Supply ability	10,000pcs/day, 1*40 HQ container can be finished in 15 days
Delivery time	According to the quantity, normally 30days after deposit
Packing details	Sterile: 1 Pcs/paper plastic bag,50 Pcs/Ctn Non-sterilize: 1 Pcs/bag, 50 Pcs/ctn Or as your requirements

### **Features of PP patient slide sheet**

- 1) Made of high quality nonwoven fabrics, non woven SMS,PP, or PP+PE 60g-120g;
- 2) With handle, easily used for transferring the patient body;
- 3) Supports up to 300 pounds. Size: 85cm W x 195cm L

### **Picture of PP heavy weight patient slide sheet**



orient Honest  
MEDICAL INTERNATIONAL

