

## Product details

### Parameter masker medis anak

- 1) Bahan: Kertas tinggi filtrasi, PP non woven + PP non woven, sangat lembut, tidak iritasi.
- 2) warna: Putih
- 3) telinga loop: loop elastis datar telinga, lebar 6 mm
- 4) ukuran: 125 × 80 mm
- 5) hidung bar: tidak ada bar hidung
- 6) BFE (bakteri efisiensi filtrasi) > 99%

### Kemasan nonwoven medis anak masker

		Ukuran karton	1 * 40 HQ kontainer dapat menahan (Pcs)
Normal Kemasan	50 Pcs/kotak, 2000 Pcs/Ctn	52 * 29 * 34 cm	2,620,000 Pcs
Satu Kemasan	1 Pcs/1 warna dicetak disegel tas plastik, 50 Pcs/box, 2000 Pcs/Ctn	52 * 30 * 42 cm	2,060,000 Pcs

### Pictures anak medis pakai masker



orient Honest  
MEDICAL INTERNATIONAL

