

All gowns are made of highest quality material, providing excellent protection and superior comfort:

Excellent material, waterproof and it's easy to use,

With thumb loop, comfortable, and easy to work, with ties behind waist, can be adjustable according to different circumference, fits your body more better.

Various size can be customized according to different requirement

Fabric material	CPE
Style	Long sleeves with thumb, with ties behind waist
Dimension sizes	S-5XL
Fabric weight	30-70g, can be as your request
Color	white, baby blue
Certificate	ISO 13485, CE, FDA
MOQ	100000 PCS
Packing	1pcs/bag, 10pcs/middle bag, 200 pcs /carton
Capacity	60000 pcs /day
Delivery time	7-30 days
Application	laboratory, clinic, food industry, dust-free workshop
Packing design	can be customized as your design



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

