

All boot cover are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Various Style can be customized according to your different requirement

| | |
|------------------------|--|
| Fabric Material | Microporous |
| Style | With elastic on ankle |
| Dimensions size | 40*50 cm can be customized according to your request |
| Fabric weight | 50-60g, could be as your request |
| Color | White color |
| Certificate | ISO13485, FDA,CE |
| MOQ: | 10000 pcs |
| Packing design | can be as your design |
| Packing | 100pcs/bag,1000 pcs/carton, or as your requirement |
| Capability | 60,000pcs per day |
| Delivery time | 7-30 days |
| Applications | laboratory, clinic, food industry, dust-free workshop,hospital |

The picture of [disposable boot cover medical manufacturer](#)







