

Product details

Indications of tie back surgical mask

- 1) Sterile face mask is for hospital, clinical, medical personnel and patients;
- 2) Non-sterile face mask can be used for dust-free workshop, food factory and so on.

Properties of protective face mask

- 1) High filtration: 3-ply pleated design offers 0.05 micron filtration at greater than 99%.
- 2) Slim aluminum wire fits comfortably over the nose, and maintain maximum filtration at the same time.
- 3) Tie and ear loop style: two pairs of polypropylene ties insure a comfortable fit;
ear loop style attaches fast and with elastic loops.

Packaging of medical mouth covers

| Specification | Normal packing | Box size | Carton size |
|---|--------------------------|------------|-------------|
| 17.5*9.5cm, 3ply | 50 pcs/box, 40 boxes/ctn | 18*10*10cm | 52*38*42cm |
| Remark: all can be special customized according to your actual requirements | | | |

Product display



The regular shipment way



FAQ

Q1: What is your output capacity for your major range of the product?

A: Face mask capacity: 1000,000 pieces / day.

Coverall, Surgical gowns, lab coat, bed sheet capacity: 10,000 pieces/day.

Caps capacity: 200,000 pieces/day. Each 1*40HQ can be finished within 1 week.

Q2: what is your MOQ of each products?

A: Face mask/caps/shoe covers/arm sleeves each MOQ: 50,000 pieces.

Clothes and sheets MOQ:5000 pieces