

lightweight comfortable [disposable bed cover for hospital /SPA](#)

We are [disposable bed cover manufacturer in China](#) .Made from high-performance isolation protective material. Widely used in hospital for emergency bed and surgical bed . [disposable bed cover](#) is used in hotels and massage bed . To prevent the spread of body fluids and blood splash content and isolate dust ,particle , alcohol, blood , bacteria and virus invading .

Material	PP,PP/PE,SMS,SMMS
Style	Flat, elastic ends, square ends, roll
Size	80x190+20cm,90x200+20cm,160x190+20cm,160x200+20cm,180x200+20cm.....
Color	white,green,blue,yellow,pink etc.
Packing	10pcs/box,10 bags/ctn or 5 pcs/bag,50 pcs/bag
Packing Design	can be as your design
MOQ	5000pcs
<b>Samples</b>	can be supplied free for your quality checking within 3 days
<b>Delivery time</b>	Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days
<b>FOB port</b>	Wuhan/Shanghai/Ningbo/Guangzhou

Picture





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL

