

La proprietà di Disposable White PP lab coat

Tutti cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:

1. Size: S-5xl or customized
2. Made of nonwoven Materiale, breathable and comfortable, clean and sanitary.
3. Convenient and lightweight , easy to carry out .

The Specification of Disposable White PP lab coat

Nome	Disposable White PP lab coat
Materiale	pp
Stile	Knitted collar and knit cuffs, with bottons and pockets
Misurare	S-5xl
Colore	blu
Imballaggio	5pcs/box,10 bags/ctn
packing design	all inner box and carton printing can be as your design
Moq	5000 pcs
Campione	can be supplied free for your quality checking within 2 days
Payment way	T/T, PayPal, Western Union, L/C,D/A,D/P,ESCROW,
Tempo di consegna	Capacity: 10,000pcs/day, 1*40HQ can be finished within 20-25 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

The picture of Disposable White PP lab coat





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

