

Handwritten header text

Handwritten text

- 1. Non : 50 pcs/40boxes/ctn 52x38x30cm
 - 1 * 20 : 500 (1,000,000)
 - 1 * 40HQ : 1150 (2,300,000 pc)
- 2. Sterile: 52x19x30cm, 1000 /ctn 50bags/ 1 /

Handwritten text

- 1. Non
- 2. Non
- 3. Non

Handwritten text

Material	disposable
Model	-2001
Structure	3ply PP + PP
Size	17.5x9.5cm(adult)
Feature	loop(latex free)
Structure	3ply
Material	
Quantity	50 /40boxes/ctn 2000pcs 1 * 20 : 500 (1,000,000) 1 * 40HQ : 1150 (2,300,000 pc)
Dimensions	52 * 38 * 30 cm
GW	GW: 7.5 kg/ctn NW: 7 kg/
MOQ	5000
Standards	CE/ISO13485/FDA

Handwritten text





orient Honest
MEDICAL INTERNATIONAL



