

## Product details

### SMS

1. SMS
2. SMS electricit
3. SMS
4. SMS
5. SMS

### SMS

	SMS
	SMS +
	S * 130 M 115 * 137 L 120 * 140 XL 125 * 150 2XL 130 * 150
	16-65gsm
	10 / 50 / CTN 100 pcs / CTN 1 / 50 / CTN 100 / CTN
	& 038; / /
	One 40HQ 15-25days finished
	80
	T / T L / C D / A D / P
	1. SMS
	2. Validity 3

OEM	1. Material <input type="checkbox"/> customers' requirements <input type="checkbox"/> 2. Customized <input type="checkbox"/> / <input type="checkbox"/> <input type="checkbox"/> 3. Customized <input type="checkbox"/>
<input type="checkbox"/>	CE <input type="checkbox"/> ISO13485 <input type="checkbox"/> FDA

## Product display



orient Honest  
MEDICAL INTERNATIONAL






orient Honest  
MEDICAL INTERNATIONAL



## Our business process



## The regular shipment way



## FAQ

**Q1. Are you a factory?**

A: Yes, we are a factory. We have been producing since 2003 and have a factory area of 2013 square meters.

**Q2. What is your price range?**

A: Our price range is USD 18 to 20 per unit.

**Q3. Do you have any certifications?**

A: Yes, we have ISO 9001, EN 15189, CE, ISO 13485, FDA, REACH, and NELSON certifications.

**Q4. Can you provide samples?**

A: Yes, we can provide samples. The lead time for samples is 3-7 days.