

## Product details

Material: Polypropylene (PP)

1. Weight: 40g
2. Material: Polypropylene
3. Material: Polypropylene
4. Size: S-XXXL
5. Quantity: 5PCS / 50PCS / 500PCS
6. Dimensions: 1 \* 20 "FCL 50000PCS
7. Quantity: 5000PCS
8. Certifications: CE, FDA, ISO13485

Material: Polypropylene (PP)

Material	PP
Weight	40g
Material	Polypropylene
Size	S-6XL
Quantity	5PCS / 50PCS / CTN
Dimensions	1 * 20 "FCL 50000PCS
Quantity	5000PCS
Certifications	dustfree, CE, FDA, ISO13485

## Product display



orient Honest  
MEDICAL INTERNATIONAL






orient Honest  
MEDICAL INTERNATIONAL



## Our business process



## The regular shipment way



## FAQ

Q1. Are you a manufacturer or a trading company?

A. We are a manufacturer since 2003. We have a factory in Shenzhen, China.

Q2. What is your minimum order quantity (MOQ)?

A. Our MOQ is 1000 units. For some products, we can accept smaller orders, but the price will be higher. The price range is USD 18 to 20 per unit.

Q3. Do you have any certifications or standards?

A. Yes, we have several certifications and standards, including ISO 9001, CE, ISO 13485, FDA, REACH, and NELSON.

Q4. Can you provide samples for testing?

A. Yes, we can provide samples for testing. The sample fee is 3-7 USD per unit, which is refundable upon a large order.