

Product details

- 1 Spunbonded SMS 20-65gsm
- 2 S - 6XL
- 3 1000/1000 / CTN

Table with 2 columns: Attribute, Value

Material	SMS
Weight	SMS40G
Color	
Size	SMS40G
Quantity	
Size	S-6XL
Quantity	1000/1000 / CTN
Material	
Weight	20-65gsm

- 1
- 2
- 3 L / C/T / T/D / P;
- 4. FOB / CNF / CIF / EXW
- 5
- 6
- 7 OEM
- 8 25-30 40HQ
- 9
- 10

Product display



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INSTRUMENTAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q1. Are you a factory?

A: Yes, we are a factory. We have been producing since 2003 and have a factory area of 2013 square meters.

Q2. What is your price?

A: Our price is very competitive, ranging from USD 18 to 20 per unit.

Q3. Do you have any ISO or EN or CE certification?

A: Yes, we have CE, ISO 13485, FDA, REACH, and NELSON certifications.

Q4. Can you deliver to my country?

A: Yes, we can deliver to your country within 3-7 days.

The happy time with customers

