

Product details

Material: PE coated. With PE coating on the inside and outside of the bag.

Weight: 100g / 50g / 50g / CTN

Material: PP + PE45G

Weight: 45gsm

Material: PP + PE45G

Weight: 45gsm

Material: PP + PE45G

Material: PP + PE45G

Material: S / M / L / XL / XXL / XXXL

Material: PP + PE45G

Material: 1PCS / 50g / 50g

Material: 80gsm

Material: CE ISO13485 FDA

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q1. Are you a manufacturer or a trading company?

A: We are a manufacturer since 2003. We have our own factory and R&D department.

Q2. What is your minimum order quantity (MOQ)?

A: Our MOQ is 1000 units. For some products, we can accept orders as low as USD 1800-2000.

Q3. Do you have any certifications or standards?

A: Yes, we are certified by CE, ISO 13485, FDA, REACH, and NELSON.

Q4. Can you provide samples for testing?

A: Yes, we can provide samples for testing. The cost of samples is 3-7 USD per unit.