

XXXXXXXXXXXXXXXXXXXX

XXXXXXXX ORIENT HONEST

XXXXXXXX OH-1001/OH1002/OH1003

XXXX CE / ISO13485 / FDA

XXXX XX

XX PPXXXX

XX XXX


XXXXXXXXXXXXXXXXXXXX

XXXX 10XXXX30XXXX

XXXX 18 "19 "21 "XX

XXXXXXXXXXXX 100/XX20/ CTN

XXXXXXXXXXXXXXXXXXXX



XX19 " 20 " 21 "XXXXXXXXXXXXXXXXXXXXXXXXXXXX

XXXX PPXXXX

XX 10XXXX16XXXX

XXXX XXXXXXXXXXXXXXXXXXXXXXX

XXXX 19 "XXXX18 "20 "21 "XX

XX XXXXXXXXXXXXXXX XX XX XX

MOQ 100000


XXXX 100/XX 2000PCS / CTN

XXXX XXX XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX XXXXXXXXXXX XXXXXXX

XXXXXXXX XXXXXXXXXXX80XXXXXXXXXXXXXXXXXXXX

XXXX T / T/L / CXXXXD / P/D / AXXXX

XXXXXXXXXXXX OEMXXXXXXXXXXXX/XXXX





orient Honest
MEDICAL INTERNATIONAL

