

**CE ISO** 00 00 00 0000 0000 00 00 00 00 0

0 0000 00

- 1) 00 00 0000 00 0000 0000 00 fe 00 99%.
- 2) 0 0000 00 3 0000, 00 0000 00000 000000 0 0 00 00 0.
- 3) 00 00 0 00 00, 00 00 0000 0 0000-000000. 00 000000 0 00 0000 00000 00 00 0 0000 0000.
- 4) 0 00 00 00 0 00: 145 × 90mm

0000 0000 00 00 0000

00	00000 00	00 00	1 * 40 HQ 000000 (Pcs) 0 00 0 00
3 00, 25 + 25 + 25g	50 pc 00 00, 40 00 00/ctn	52 * 31 * 34 cm	2440 Pcs
0: 0000 00 00 0 0000 0000 00000 0 00 0000 000000 00 0000			

00 00 0000 00 00000 00 0000

- 1) 00000 0
- 2) 00 0000 00 00
- 3) 000000 00 0 00, 0 0000 0000 0 0 00 00 00000 000000
- 4) 00000000 00 0 00, 0 0000 0000 0 0 00 00 000000 000000.
- 5) 0000 00 0 0000 000000.

00 00000 00



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL

