

□□ □□ □ □ □□ □□ □ □□ □□ □□

□□□□ □□

- 1. Non □ □: 50pcs/□□, 40boxes/ctn, 52x38x30cm.  
 1 \* 20 □□ □ □: 500□ □□ (1000,000 pc),  
 1 \* 40HQ □□: 1150 □□□□ (2300,000 pc) □ □ □□.
- 2. Sterile: 1pcs/□□, 50bags/□□, 1000pcs/ctn □□ □□: 52x19x30cm.

□□

- 1. □ □ □□□ □□ □□ □□□ □ □□ □ □ □□ □□.
- 2. Non-□ □ □□□ □□-□□ □□□, □□, □□□, □□□ □ □□□□.
- 3. □□□ □□□ □□ □□ □□□ □□ □ □□ □□□□ □□□ □□ □□,

□□

|       |  |
|-------|--|
| □□ □□ | □□ disposable □ □ □□ □□ □ □□ □□ □□ □□ □□   |
| □□ □□ | □-2001   |
| □□    | 3ply, PP □ □ + PP □ □ □□ □□  |
| □□    | 17.5x9.5cm(adult)  |
| □□□   | □□ □□ □-loop(latex free)   |
| □□□   | 3ply   |
| □     | □□□  |
| □□    | 50pcs/□□, 40boxes/ctn, 2000pcs/ctn.<br>1 * 20 □□ □ □: 500□ □□ (1000,000 pc),<br>1 * 40HQ □□: 1150 □□□□ (2300,000 pc) □ □ □□. |
| □□ □□ | 52 * 38 * 30 cm  |
| GW    | GW: 7.5 k g/ctn, NW: 7kg/ctn.  |
| MOQ   | 5000PCS  |
| □□□   | CE / ISO13485 / □□ / FDA   |

□□





orient Honest  
MEDICAL INTERNATIONAL



