

3 ply surgical face mask , anti-virus breathing mask , single use dental face mask

Parameters of 3 ply surgical face mask

- 1) Material: PP non woven + high filtration paper + PP non woven
- 2) Ear loop: tie on or elastic ear loop
- 3) Nose bar: plastic or metal, flexible, fitting fixed
- 4) Size: 17.5*9.5 cm, one size fit all.
- 5) Color: White/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/back

Packaging of breathing mask

| Specification | Normal packing | Carton size | 1*40 HQ container can hold (Pcs) |
|---|-----------------------------|-------------|----------------------------------|
| 17.5*9.5cm,3ply, 18+20+20G | 50 pcs/box, 40 boxes/ctn | 52*19*42cm | 3,280,000 Pcs |
| 17.5*9.5cm,3ply, 25+25+25G | 50 pcs/box, 40 boxes/ctn | 51*37*30cm | 2,400,000 Pcs |
| Remark: all can be special customized according to your actual requirements | | | |

Our advantage

- 1) We are manufacturer for dental face mask with 12 years history, have large advantage in the price and technology;
- 2) Strong production power:1*40 HQ container can be finished in 15 days;
- 3) Promising sales team: service in 24 hours;
- 4) Sample is free, sample can be finished in 3 days;
- 5) High quality insurance: ISO9001, ISO13485, CE certificated, FDA Registered.

Picture of 3 ply surgical face mask

